

We provide the best

for your laboratory needs

along this path









We offer the best solution in pathology for better diagnosis

•Tissue Processor •Paraffin Dispenser •Cold Plate •Paraffin Trimmer •Rotary Microtome •Tissue Float •Hot Plate













Tissue Processor

Pooyan Mk2250

- ·Colorful touchscreen in persian and english
- •SMS module to send errors & monitor specimens situation
- Rechargable battery to protect specimens during a power failure from tissues drying
- Stirring system for better process
- •Gentle speciemen processing with high safety for tissue
- ·User friendly interface design
- ·Ability to start the program from any of the station
- •Carousel type with 12 stations of 1.8 liters each;10 reagent stations, 2 wax bath
- •Two stainless steel tissue basket, capacity: approx.100 cassettes
- •Easy editing of programs, even during a processing run
- Delayed start function up to 70 days





Paraffin Dispenser

Pooyan Mk3335

- Adjustable timer to turn the device on and off
- LED lamp for lighting work area
- ·Ability to adjust paraffin reservoir & work surface temperature
- User friendly control panel for easy setup of parameters
- Activate the paraffin flow via manual switch or foot pedal
- Maximum efficiency and operating comfort
- Forceps retaining and heating holes suitable for left or right hand operation
- ·Large paraffin reservoir with filter
- Tray for the collection of overflowing paraffin



Semi Automated Rotary Microtome

Pooyan Mk3130

- •Section thickness setting range from 0.5 100 μm
- •Trimming thickness setting range from 5 100 μm
- •Colorful display to show section counter, trim & section thickness
- •Retraction system to extend blade life & serial cutting
- •Blade holder clamp with angular, longitudinal and transverse movements
- ·Handwheel locking system & finger guard
- Speciemen oriention in X,Y,Z
- ·Smooth-runnig handwheel
- Ergonomic, compact & user friendly design
- Spacious and easy to clean section waste tray
- •Automatic movement of the block clamp with warning alarm at the beginning and end of the course













Microtome blade holder

Standard & universal clamp

Cold Plate

Pooyan Mk4615

- Large cooling surface for up to 60 cassettes or molds
- Steady temperature distribution over the entire surface
- •Selectable temperature between 0 and -15 °C
- Adjustable timer to turn the device on and off
- •Low power consumption (LPC) compressor
- •User friendly control panel for easy setup of parameters
- •Removable stainless steel drainage tray located under the plate
- •fully efficiency at the same time as beauty and minimum dimension



POOYAN TEB KHADEM

best solotion in pathology

Tissue Float

Pooyan Mk3420

- Removable glass dish for easy cleaning
- ·LED light array for optimum specimen viewing
- Precise temperature regulation to prevent tissue imperfections
- Heat protector to prevent overheating
- Adjustable timer to turn the device on and off
- •Selectable temperature between 30 and 65 °C



Hot Plate

Pooyan Mk4510

- Suitable surface size to dry slides
- ·Adjustable timer to turn the device on and off
- •User friendly control panel for easy setup of parameters
- •Digital control panel to set and monitor the temperature of the heating surface
- Heat protector to prevent overheating
- •Selectable temperature between 30 and 80 °C



Paraffin Trimmer

Pooyan Mk1711

- •Paraffin trimmer melt excess paraffin on outer of Tissue cassette and user can reuse paraffin to mold tissue
- •The surface of device has standard grooves for tissue cassettes and appropriate temperature which speed up your work
- •The risk of getting clean text or barcodes written on tissue cassettes is zero



Production and quality control procedures certified under EN ISO 13485
IEC 61010-1 & IEC 60601-2 ensure highest quality and reliability



Processing / embedding cassette

Pooyan Teb Khadem Co.

Add: No99, Part4, Toos industrial Park

Mashhad,iran

TEL: +98 51 3541 8561

Mobile: +98 915 110 2589 , +98 915 122 2488

Email: mkhadem.co@gmail.com Web: www.pooyantebkhadem.com







pooyan_teb_khadem